

Art. Nr.: WM-STDA-03-EN
Indication: Patients with diabetes mellitus
Location/year: Movement Analysis Lab, University Hospital Münster; Germany/2009
Surgeons: Arne Nagel, Dieter Rosenbaum
Tested shoes: VACODiaped, VACODiaped plus¹, DARCO light, "Health Shoe"



Contents lists available at ScienceDirect

Gait & Posture

journal homepage: www.elsevier.com/locate/gaitpost



Vacuum cushioned removable cast walkers reduce foot loading in patients with diabetes mellitus

Arne Nagel ^{a,b}, Dieter Rosenbaum ^{a,*}

^a Movement Analysis Lab, Orthopaedic Department, University Hospital Münster, Germany

^b Motion Analysis Lab Walstedde – Centre for Motion Analysis and Therapy, Drensteinfurt-Walstedde, Germany

Summary:

The purpose of this study was to investigate the pressure-relieving effects of two vacuum orthoses (VACOped Diabetic & VACOpedes Diabetic) in patients with diabetes mellitus. Twenty patients: 13 females & 7 males, average age 56,4 +/- 15,7 years. Disease duration was 12,8 +/- 12,1 years. They presented with plantar callosities, but no ulceration.

Method & Materials:

Plantar pressure distribution was measured with sensor insoles during walking in two different VACO-orthoses, a postoperative shoe and a common "Health Shoe".

Results:

Significant differences in plantar pressure distribution were found between the four walking conditions. The contact area increased in the midfoot with the vacuum orthosis (up to 45,1%).

→ **Maximum force & peak pressures showed a significant decrease under the rearfoot and forefoot and increased in the midfoot area during walking with both vacuum orthoses (VACOped Diabetic & VACOpedes Diabetic)!**

→ **The VACODiaped revealed equal pressure relief under the forefoot and significantly lower rearfoot pressures compared to the post-operative shoe DARCO light.**

The vacuum orthosis demonstrated a comparable pressure relieving efficacy under the forefoot to post-operative shoes.

Conclusion:

→ **For patients using a VACOped Diabetic or the VACOpedes Diabetic the maximum force & peak pressures showed a significant decrease under the rearfoot and forefoot and they significantly benefited re-distribution of plantar pressure and the roll-over process.**

Full Text Link:

<https://pubmed.ncbi.nlm.nih.gov/19321342/>

¹Old versions are called VACODiaped and VACODiaped plus

New Versions are called VACOped Diabetic and VACOpedes Diabetic