

Order direct at
Prescribe through **patientchoice**
Free on-site stock available


VACO[®]cast DIABETIC

NICE and IWGDF Offloading Guidelines recommend a non-removable knee-high walker as 1st line treatment to offload DFUs*¹

**Drug
Tariff
listed**

VACO12

Vacuum cushion
VACO12 Technology provides TCC-like offloading and shear stress protection, which can be repeatedly adjusted to the individual.



Please scan to view animation video on youtube

Skin friendly velour liner
Additional liner included. Wash one liner while wearing the other.



Optional OPED Locking Mechanism
for a non-removable device to ensure compliance.



Removable toe protector



Removable rocker sole
Wear one sole for inside and one for outside the house.



Calf extension
adaptable for extra large calf circumference

VACOcast Diabetic is in Part IX A of Drug Tariff
VACOcast Diabetic Small
Order code: VPCD-01-1-S
PIP code: 413-0522
VACOcast Diabetic Medium
Order code: VPCD-01-1-M
PIP code: 413-0530
VACOcast Diabetic Large
Order code: VPCD-01-1-L
PIP code: 413-0514
Strap Lock
Order code: VERS-16-4-41
PIP code: 418-9916



Please scan to view application video on youtube

For appointments please call **01380 722177** and for more information or application-videos please visit www.offloading-device.com.

*1 <https://www.nice.org.uk/guidance/ng19/chapter/Recommendations#diabetic-foot-ulcer>
*2 IWGDF Offloading Guideline 2019 <https://iwgdfguidelines.org/offloading-guideline>

A complete solution with no hidden costs

Each device contains:

- ready to use VACOcast Diabetic boot
- replacement liner – wash one whilst wearing the other
- replacement rocker soles – one to wear outside the house and one to wear inside for hygiene purposes
- replacement toe cover
- set of tamper seals – to identify if the patient has removed the device
- air extractor bulb

