


# New Consensus Document recommends VACOcast Diabetic as comparable offloading to TCC



## Load redistribution and load share devices: Off-the-shelf walkers (also known as removable cast walkers/instant TCC)

Vacuum cast (VACOcast<sup>®</sup> Diabetic)



- Suitable for patients with neuropathy, arterial diseases, highly exuding wounds and infected wounds.
- Not suitable for patients with leg deformities or chronic Charcot foot.\*1)
- Non-removable by the patient. The device can be locked and unlocked with a key by authorised personnel, such as community health professionals for wound dressing change.
- Available in three sizes and comes as a complete set, including replacement liner and sole.
- Provides total contact offloading comparable to TCC (Götz et al, 2017) with no pressure applied.

[Link to Article: ESS20\\_CD\\_TCC\\_DFJ\\_web.pdf \(essity.co.uk\)](https://essity.co.uk/ESS20_CD_TCC_DFJ_web.pdf)

The Diabetic Foot Journal (April 2021)

\*1) If Patients leg fits into the boot, VACOcast Diabetic can be applied also with leg deformities or Chronic Charcot. The same applies for any other pre-fabricated boot